

## Pressure Care

### PCS701 Foam overlay

- Convoluted foam overlays redistribute the body weight to aid relaxation and to relieve pressure



### PCS702 Spenco bed pad

- Hollow core fibres and a breathable cotton/polyester flame retardant cover
- The silicone coating lubricates each fibre to reduce friction while preventing moisture penetration or retention
- One side is breathable, opposite side is waterproof. Dries in minutes



### PCS703 Bennu alternating pressure overlay



## Pressure Care

### PCS706

#### Softform Premier Active mattress

- The unique design features an alternating air inset beneath the castellated foam insert of the mattress
- A pump will assess individual patient's weights and adjust levels of inflation accordingly  
Air cells alternate on a ten minute cycle beneath the foam inert, delivering a very low impact at the point of change



### PCS708

#### Sheepskins

- Medical grade sheepskins
- Help to reduce the risk of pressure sores and can be used in beds, wheelchairs or chairs
- High temperature ones available in green



### PCS707

#### JFC Foam Mattress

- The Bayer viscoelastic foam contours to the body and relieves pressure

